

Product Datasheet

CDK5 (pY15) antibody (orb1422844)

Description	Rabbit polyclonal antibody to CDK5 (pY15).
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q00535
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), IF: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500), E: WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.