

Product Datasheet

STAT6 (pY641) antibody (orb1423034)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to STAT6 (pY641). |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Porcine, Rabbit, Rat, Virus |
| Conjugation | Unconjugated |
| Tested Applications | IH, IP, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | P42226 |
| Dilution Range | WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100) |
| Expiration Date | 12 months from date of receipt. |