

Product Datasheet

PDPK1 (pS241) antibody (orb1423434)

Description	Rabbit polyclonal antibody to PDPK1 (pS241).
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Virus
Conjugation	Unconjugated
Tested Applications	IH, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O15530
Dilution Range	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200), IP (1/10 - 1/100)
Expiration Date	12 months from date of receipt.