

## Product Datasheet

### CDK2 (pT160) antibody (orb1423581)

<b>Description</b>	Rabbit polyclonal antibody to CDK2 (pT160).
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Porcine, Rabbit, Rat, Virus, Zebrafish
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IH, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P24941</b>
<b>Dilution Range</b>	WB: WB (1/500 - 1/1000), IH (1/100 - 1/200), IH: WB (1/500 - 1/1000), IH (1/100 - 1/200)
<b>Expiration Date</b>	12 months from date of receipt.