



Product Datasheet

SOCS1 antibody (orb1427740)


Description	Rabbit polyclonal antibody to SOCS1.
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	SOCS1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23551
Uniprot ID	O15524
NCBI	NP_003736.1
Dilution Range	IF: 1:10~50, WB: 1:2000, WB: 1:1000, IHC-P: 1:10~50, FC: 1:10~50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of S...



Western blot: All lanes: Anti-SOCS1 Antibody...



Western blot analysis of SOCS1 Antibody ...