



Product Datasheet

ALOX15 antibody (orb1428334)


Description	Rabbit polyclonal antibody to ALOX15.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ALOX15
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74804
Uniprot ID	P16050
NCBI	NP_001131.3
Dilution Range	IF: 1:10~50, WB: 1:2000, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of A...



Western blot:All lanes:Anti-ALOX15 Antib...



Western blot analysis of HeLa(lane 1), A...