



Product Datasheet

LIF antibody (orb1429401)


Description	Rabbit polyclonal antibody to LIF.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	LIF
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22008
Uniprot ID	P15018
NCBI	NP_002300.1 , NP_001244064.1
Dilution Range	WB: 1:1000, WB: 1:1000, IHC-P-Leica: 1:500, FC: 1:25
Expiration Date	12 months from date of receipt.



Western blot:Anti-LIF Antibody (Center) ...



Western blot:All lanes:Anti-LIF Antibody...



Immunohistochemical analysis of human ki...