

Product Datasheet

RAGE (AGER) antibody (orb1431608)

Description	Rabbit polyclonal antibody to RAGE (AGER).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	AGER
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42803
Uniprot ID	Q15109
NCBI	NP_001193895.1 , NP_001193863.1 , NP_001193869.1 , NP_001193858.1 , NP_001193883.1 , NP_001193865.1 , NP_001193861.1 , NP_751947.1 , NP_001127.1
Dilution Range	WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.