

Product Datasheet

ANKHD1 antibody (orb1435632)

Description	Rabbit polyclonal antibody to ANKHD1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ANKHD1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	269458
Uniprot ID	Q81WZ3
NCBI	NP_065741.3 , NP_060448.1 , NP_001183959.1 , NP_078944.2 , NP_060217.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.