



Product Datasheet

HSPA5 antibody (orb1437589)


Description	Rabbit polyclonal antibody to HSPA5.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	HSPA5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72333
Uniprot ID	P11021
NCBI	NP_005338.1
Dilution Range	IF: 1:10~50, WB: 1:1000, WB: 1:1000, IHC-P: 1:50~100, IHC-P: 1:50~100, FC: 1:10~50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



The anti-HSPA5 Pab is used in Western bl...



Western blot analysis of anti-HSPA5 Pab ...