

Product Datasheet

INA antibody (orb1440596)

Description	Mouse monoclonal antibody to INA.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	INA
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	55391
Uniprot ID	Q16352
NCBI	NP_116116.1
Dilution Range	IF: 1:10~50, WB: 1:125, WB: 1:1000, IHC-P: 1:10~50
Expiration Date	12 months from date of receipt.