

## Product Datasheet

### p38 MAPK (Phospho-Y182) Blocking peptide (orb1441116)

|                          |                                                                                                                                      |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | p38 MAPK (Phospho-Y182) Blocking peptide                                                                                             |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| <b>Note</b>              | For research use only                                                                                                                |
| <b>Application notes</b> | Synthetic peptide                                                                                                                    |
| <b>Protein Sequence</b>  | The synthetic peptide sequence is selected from aa 175-186 of HUMAN MAPK11                                                           |
| <b>Uniprot ID</b>        | <b>Q15759</b>                                                                                                                        |
| <b>Expiration Date</b>   | 12 months from date of receipt.                                                                                                      |