

Product Datasheet

ARHGAP35 (Phospho-Y1105) Blocking peptide (orb1441119)

Description	ARHGAP35 (Phospho-Y1105) Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 1090-1115 of HUMAN ARHGAP35
Uniprot ID	Q9NRY4
Expiration Date	12 months from date of receipt.