

## Product Datasheet

# CCND1(Phospho-T288) Blocking peptide (orb1441138)

<b>Description</b>	CCND1(Phospho-T288) Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence is selected from aa 282-294 of HUMAN CCND1
<b>Uniprot ID</b>	<b>P24385</b>
<b>Expiration Date</b>	12 months from date of receipt.