

Product Datasheet

CCND3 (Phospho-T283) Blocking peptide (orb1441152)

Description	CCND3 (Phospho-T283) Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 277-288 of HUMAN CCND3
Uniprot ID	P30281
Expiration Date	12 months from date of receipt.