

## Product Datasheet

# YAP (Phospho-S127) Blocking peptide (orb1441242)

<b>Description</b>	YAP (Phospho-S127) Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence is selected from aa 119-133 of HUMAN YAP1
<b>Uniprot ID</b>	<b>P46937</b>
<b>Expiration Date</b>	12 months from date of receipt.