

## Product Datasheet

### ZDHHC17 Blocking peptide (orb1441502)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | ZDHHC17 Blocking peptide   |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Synthetic peptide  |
| <b>Protein Sequence</b>  | The synthetic peptide sequence is selected from aa 3-17 of HUMAN ZDHHC17   |
| <b>Uniprot ID</b>        | <b>Q8IUH5</b>  |
| <b>Expiration Date</b>   | 12 months from date of receipt.  |