

Product Datasheet

SSX1 Blocking peptide (orb1441545)

| | |
|--------------------------|--|
| Description | SSX1 Blocking peptide |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Note | For research use only |
| Application notes | Synthetic peptide |
| Protein Sequence | The synthetic peptide sequence is selected from aa 117-128 of HUMAN SSX1 |
| Uniprot ID | Q16384 |
| Expiration Date | 12 months from date of receipt. |