

Product Datasheet

Ehmt2 Blocking peptide (orb1442158)

Description	Ehmt2 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 227-239 of HUMAN Ehmt2
Uniprot ID	Q9Z148
Expiration Date	12 months from date of receipt.