

Product Datasheet

(DANRE) gabarapl2 Blocking peptide (orb1442226)

Description	(DANRE) gabarapl2 Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 15-28 of HUMAN Q6TH05
Uniprot ID	Q6TH05
Expiration Date	12 months from date of receipt.