

## Product Datasheet

### ECE2 Blocking peptide (orb1442942)

<b>Description</b>	ECE2 Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence is selected from aa 566-579 of Human ECE2
<b>Uniprot ID</b>	<b>O60344</b>
<b>Expiration Date</b>	12 months from date of receipt.