

Product Datasheet

GABARAPL2 (GATE16) Blocking peptide (P71) (orb1443550)

Description	GABARAPL2 (GATE16) Blocking peptide (P71)
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 71~82 of HUMAN GABARAPL2
Uniprot ID	P60520
Expiration Date	12 months from date of receipt.