

## Product Datasheet

### JHDM1a Blocking peptide (orb1443599)

<b>Description</b>	JHDM1a Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence is selected from aa 882-896 of HUMAN KDM2A
<b>Uniprot ID</b>	<a href="#">Q9Y2K7</a>
<b>Expiration Date</b>	12 months from date of receipt.