

## Product Datasheet

### NKX3-1 Antibody Blocking peptide (orb1444637)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | NKX3-1 Antibody Blocking peptide  |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.  |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Synthetic peptide   |
| <b>Protein Sequence</b>  | The synthetic peptide sequence used to generate the antibody orb1928164 was selected from the Center region of human NKX3-1. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| <b>Uniprot ID</b>        | <a href="#">Q99801</a>  |
| <b>Expiration Date</b>   | 12 months from date of receipt.   |