

Product Datasheet

PKN beta Antibody Blocking peptide (orb1445673)

Description	PKN beta Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1928751 was selected from the C-term region of human PKN beta . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q6P5Z2
Expiration Date	12 months from date of receipt.