

Product Datasheet

CYP21A2 Antibody Blocking peptide (orb1445761)

| | |
|--------------------------|--|
| Description | CYP21A2 Antibody Blocking peptide |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Note | For research use only |
| Application notes | Synthetic peptide |
| Protein Sequence | The synthetic peptide sequence used to generate the antibody orb1928798 was selected from the N-term region of human CYP21A2. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| Uniprot ID | P08686 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)