



Product Datasheet

STAT3 (Ab-705) Antibody (orb14458)


| | |
|----------------------------|---|
| Description | STAT3 (Ab-705) Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P40763 |
| Expiration Date | 12 months from date of receipt. |



Immunohistochemical analysis of formalin...



Western blot analysis of STAT3(Ab-705) a...



Confocal immunofluorescent analysis of H...