



Product Datasheet

STAT5A (Ab-780) Antibody (orb14462)

Description	STAT5A (Ab-780) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P42229
Expiration Date	12 months from date of receipt.



Western blot analysis of STAT5A(Ab-780) ...



Immunohistochemical analysis of formalin...