

Product Datasheet

JLP (SPAG9) Antibody Blocking peptide (orb1446757)

Description JLP (SPAG9) Antibody Blocking peptide

Form/Appearance Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Note For research use only

Application notes Synthetic peptide

Protein Sequence The synthetic peptide sequence used to generate the antibody orb1929391 was selected from the N-term region of human SPAG9. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

Uniprot ID [O60271](#)

Expiration Date 12 months from date of receipt.