



Product Datasheet

PTEN (Ab-380/382/383) Antibody (orb14468)

Description	PTEN(Ab-380/382/383) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P60484
Expiration Date	12 months from date of receipt.



Western blot analysis of PTEN(Ab-380/382/383) antibody.



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue.