

Product Datasheet

ACVR2A Antibody Blocking peptide (orb1447080)

Description	ACVR2A Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1929553 was selected from the N-term region of human ACVR2A. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	P27037
Expiration Date	12 months from date of receipt.