



Product Datasheet

Estrogen Receptor- α (Ab18) Antibody (orb14479)


Description	Estrogen Receptor- α (Ab18) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P03372
Expiration Date	12 months from date of receipt.



Western blot analysis of Estrogen Recept...



Confocal immunofluorescent analysis of M...



Immunohistochemical analysis of formalin...