

Product Datasheet

Pael-R (GPR37) Antibody Blocking peptide (orb1447949)

| | |
|--------------------------|--|
| Description | Pael-R (GPR37) Antibody Blocking peptide |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Note | For research use only |
| Application notes | Synthetic peptide |
| Protein Sequence | The synthetic peptide sequence used to generate the antibody orb1930047 was selected from the GPR37 region of human Pael-R (GPR37). A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| Uniprot ID | O15354 |
| Expiration Date | 12 months from date of receipt. |