

Product Datasheet

C-rel Antibody(NFkB) Antibody(D504) Blocking peptide (orb1448024)

Description	C-rel Antibody(NFkB) Antibody(D504) Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence was selected from the NFkB region of human C-rel (NFkB) (C-term D504). A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q04864
Expiration Date	12 months from date of receipt.