

Product Datasheet

Kallikrein 7 Antibody Blocking peptide (orb1448053)

Description	Kallikrein 7 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1930085 was selected from the Center region of human KLK7. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	P49862
Expiration Date	12 months from date of receipt.