

Product Datasheet

ACE2 (NCOVID/SARS Receptor) Antibody Blocking peptide (orb1448375)

Description	ACE2 (NCOVID/SARS Receptor) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence was selected from the SARS region of human ACE2 (SARS Receptor) . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q9BYF1
Expiration Date	12 months from date of receipt.