



Product Datasheet

EGFR (Ab070) Antibody (orb14485)


Description	EGFR (Ab070) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of M...



Western blot analysis of EGFR(Ab-1070) a...



Immunohistochemical analysis of formalin...