

## Product Datasheet

# ULK1 (Phospho-S556) Antibody Blocking peptide (orb1449621)

<b>Description</b>	ULK1 (Phospho-S556) Antibody Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">O75385</a>
<b>Expiration Date</b>	12 months from date of receipt.