

## Product Datasheet

# TSC1 (Phospho-S555) Antibody Blocking peptide (orb1449626)

<b>Description</b>	TSC1 (Phospho-S555) Antibody Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9EP53</b>
<b>Expiration Date</b>	12 months from date of receipt.