

## Product Datasheet

# ERBB2 (Phospho-Y1223) Antibody Blocking peptide (orb1449682)

<b>Description</b>	ERBB2 (Phospho-Y1223) Antibody Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">P70424</a>
<b>Expiration Date</b>	12 months from date of receipt.