



Product Datasheet

p53 (Ab5) Antibody (orb14497)


Description	p53(Ab-15) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P04637
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



Immunohistochemical analysis of formalin...



Western blot analysis of p53(Ab-15) anti...