

## Product Datasheet

# BAD (Phospho-T80) Antibody Blocking peptide (orb1449700)

<b>Description</b>	BAD (Phospho-T80) Antibody Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">Q92934</a>
<b>Expiration Date</b>	12 months from date of receipt.