

Product Datasheet

E2F1(Phospho-S364) Antibody Blocking peptide (orb1449930)

| | |
|--------------------------|--|
| Description | E2F1(Phospho-S364) Antibody Blocking peptide |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Note | For research use only |
| Application notes | Synthetic peptide |
| Protein Sequence | The synthetic peptide sequence was selected from the region of human Phospho-E2F1-S364. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| Uniprot ID | Q01094 |
| Expiration Date | 12 months from date of receipt. |