

## Product Datasheet

# EGFR (Phospho-Y1092) Antibody Blocking peptide (orb1450052)

<b>Description</b>	EGFR (Phospho-Y1092) Antibody Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence used to generate the antibody orb1931172 was selected from the region of human Phospho-EGFR-pY1092. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
<b>Uniprot ID</b>	<b>P00533</b>
<b>Expiration Date</b>	12 months from date of receipt.