

Product Datasheet

MAP3K5(Phospho-S1033) Antibody Blocking peptide (orb1450123)

Description	MAP3K5(Phospho-S1033) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence was selected from the region of human Phospho-MAP3K5-S1033. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q99683
Expiration Date	12 months from date of receipt.