

## Product Datasheet

# Smc1(Phospho-S957) Antibody Blocking peptide (orb1450186)

<b>Description</b>	Smc1(Phospho-S957) Antibody Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence was selected from the 950-964 region of human Phospho-Smc1-S957. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
<b>Uniprot ID</b>	<b>Q14683</b>
<b>Expiration Date</b>	12 months from date of receipt.