

Product Datasheet

Rb (Phospho-S811) Antibody Blocking peptide (orb1450197)

Description	Rb (Phospho-S811) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1931227 was selected from the 804-818 CR>region of human PhosphoRb-S811. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	P06400
Expiration Date	12 months from date of receipt.