

Product Datasheet

cJun (Phospho-S63) Antibody Blocking peptide (orb1450309)

Description	cJun (Phospho-S63) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1931278 was selected from the 56-67 CR> region of human Phospho-cJun-S63. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	P05412
Expiration Date	12 months from date of receipt.