

Product Datasheet

RCE1 (W57) Antibody Blocking peptide (orb1450893)

Description RCE1 (W57) Antibody Blocking peptide

Form/Appearance Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Note For research use only

Application notes Synthetic peptide

Protein Sequence The synthetic peptide sequence used to generate the antibody orb1931614 was selected from the N-term region of human RCE1 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

Uniprot ID [Q9Y256](#)

Expiration Date 12 months from date of receipt.