

Product Datasheet

Ubiquilin1 (PLIC1) Antibody Blocking peptide (orb1450989)

| | |
|--------------------------|---|
| Description | Ubiquilin1 (PLIC1) Antibody Blocking peptide |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Note | For research use only |
| Application notes | Synthetic peptide |
| Protein Sequence | The synthetic peptide sequence used to generate the antibody orb1931659 was selected from the N-term region of human PLIC1 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| Uniprot ID | Q9UMX0 |
| Expiration Date | 12 months from date of receipt. |